

Product datasheet

Goat Anti-Guinea pig IgG H&L (Alkaline Phosphatase) preadsorbed ab7140

4 References

概述

产品名称	山羊抗豚鼠IgG H&L (Alkaline Phosphatase)预吸附二抗
宿主	Goat
靶标种属	Guinea pig
经测试应用	适用于: ICC/IF, Dot blot, ELISA, IHC-P, IHC-Fr, Immunomicroscopy, WB
最小交叉反应	Chicken, Cow, Goat, Hamster, Horse, Human, Mouse, Rabbit, Rat, Sheep more details
免疫原	guinea pig IgG whole molecule
偶联物	Alkaline Phosphatase

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Constituents: 0.88% Sodium chloride, 50% Glycerol (glycerin, glycerine), 0.0014% Zinc chloride, 0.0095% Magnesium chloride, 0.79% Tris HCl
纯度	Affinity purified
纯化说明	This product was prepared from monospecific antiserum by immunoaffinity chromatography using guinea pig IgG coupled to agarose beads.
共轭说明	Alkaline Phosphatase (Calf Intestine) (Molecular Weight 140,000 daltons)
克隆	多克隆

应用

The Abpromise guarantee

Abpromise™承诺保证使用ab7140于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ICC/IF		
Dot blot		
ELISA		
IHC-P		
IHC-Fr		
Immunomicroscopy		
WB		

应用说明

suitable for immunoblotting (western or dot blot)
ELISA and immunohistochemistry as well as other phosphatase-antibody based enzymatic assays requiring lot-to-lot consistency

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors