

Goat Anti-Human IgG (Fab')₂ (Biotin) preadsorbed ab64666

2 References

概述

产品名称	山羊抗人IgG (Fab') ₂ (Biotin)预吸附二抗
宿主	Goat
靶标种属	Human
特异性	ab64666 is specific for Fab and F(ab') ₂ portion of human IgG by IEP. It reacts with other human immunoglobulins through common light chain reactivity. No reactivity with the Fc portion of human IgG or non immunoglobulin serum proteins was detected.
经测试应用	适用于: WB, ELISA, IHC-Fr, Flow Cyt
最小交叉反应	Cow, Horse, Mouse more details
免疫原	Human IgG F(ab') ₂
偶联物	Biotin

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.05% Sodium azide Constituent: PBS
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab64666于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		1/1000 - 1/10000.
ELISA		1/1000 - 1/10000.
IHC-Fr		1/200 - 1/1000.
Flow Cyt		1/200 - 1/1000.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors