

# Goat Anti-Human IgG Fc (HRP) ab97225

## 48 References

### 概述

产品名称	山羊抗人IgG Fc (HRP)
宿主	Goat
靶标种属	Human
特异性	By immunoelectrophoresis and ELISA this antibody reacts specifically with Human IgG. Cross reactivity with other IgM, IgA, and light chains is less than 1%. This antibody may cross react with IgG from other species.
经测试应用	<b>适用于:</b> ICC, IHC-P, ELISA, WB
免疫原	Other Immunogen Type corresponding to Human IgG.
偶联物	HRP

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 6.80 Constituents: 0.2% BSA, 99% PBS, 0.05% CMIT/MIT based preservative
纯度	Immunogen affinity purified
纯化说明	Antiserum was solid phase adsorbed to ensure class specificity. This antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to Horse Radish Peroxidase (HRP).
共轭说明	Molar enzyme/ antibody protein ratio is 4:1
克隆	多克隆
同种型	IgG

### 应用

#### The Abpromise guarantee

**Abpromise™** 承诺保证使用 ab97225 于以下的经测试应用

“应用说明”部分下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ICC		Use at an assay dependent dilution.
IHC-P		1/200 - 1/2000.
ELISA		1/10000 - 1/100000. (Primary)
WB		1/2000 - 1/25000. Colormetric: 1/2000 - 1/20000; Chemiluminescent: 1/5000 - 1/25000

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors