

# Goat Anti-Human IgG H&L (Alkaline Phosphatase) preadsorbed ab97167

### 概述

产品名称	山羊抗人IgG H&L (Alkaline Phosphatase)预吸附二抗
宿主	Goat
靶标种属	Human
特异性	By immunoelectrophoresis and ELISA this antibody reacts specifically with Human IgG and with light chains common to other Human immunoglobulins. No antibody was detected against non immunoglobulin serum proteins. Reduced cross-reactivity to bovine, chicken, goat, horse, mouse, pig, rabbit and rat was detected.
经测试应用	适用于: ICC, IHC-P, ELISA, WB
最小交叉反应	Chicken, Cow, Horse, Mouse, Pig, Rabbit, Rat
偶联物	Alkaline Phosphatase

[more details](#)

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 6.8 Preservative: 0.09% Sodium azide Constituents: 0.58% Sodium chloride, 0.0095% Magnesium chloride, 0.00136% Zinc chloride, 1.19% HEPES, 0.2% BSA
纯度	Immunogen affinity purified
纯化说明	Antiserum was cross adsorbed using bovine, chicken, horse, mouse, pig, rabbit and rat immunosorbents to remove cross reactive antibodies. This antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to Alkaline Phosphatase (AP).
克隆	多克隆
同种型	IgG

## 应用

### The Abpromise guarantee **Abpromise™**承诺保证使用ab97167于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ICC		Use at an assay dependent dilution.
IHC-P		1/100 - 1/1000.
ELISA		1/2500 - 1/25000. ELISA-Primary
WB		1/2500 - 1/25000. WB-Chemiluminescent

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors