

Goat Anti-Mouse IgG gamma chain (beta-galactosidase) preadsorbed ab130781

概述

产品名称	山羊抗小鼠IgG gamma chain (beta-galactosidase)预吸附二抗
宿主	Goat
靶标种属	Mouse
特异性	ab130781 reacts with mouse IgG1, IgG2a, IgG2b and IgG3. Goat Anti-Mouse IgG gamma chain (beta-galactosidase) preadsorbed exhibits minimal cross reactivity with Human immunoglobulins.
经测试应用	适用于: ELISA, WB, ICC/IF
最小交叉反应	Human more details
免疫原	Full length protein corresponding to Mouse IgG. Mouse IgG paraproteins.
偶联物	beta-galactosidase

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C. Store In the Dark.
存储溶液	pH: 7.2 Preservative: 0.1% Sodium azide Constituents: 50% Glycerol (glycerin, glycerine), 49% PBS
纯度	IgG fraction
纯化说明	ab130781 was purified by affinity chromatography on pooled mouse IgG covalently linked to agarose.
克隆	多克隆
同种型	IgG
常规说明	ab130781 is pre-adsorbed against Mouse IgM and IgA, pooled Human sera and purified Human paraproteins.

应用

The Abpromise guarantee

Abpromise™ 承诺保证使用ab130781于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ELISA		1/500.
WB		Use at an assay dependent dilution.
ICC/IF		Use at an assay dependent dilution.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors