

Goat Anti-Mouse IgG H&L ab6708

56 References 2 图像

概述

产品名称	山羊抗小鼠IgG H&L
宿主	Goat
靶标种属	Mouse
经测试应用	适用于: ELISA, IP, Immunodiffusion, Conjugation, ICC/IF, IHC-P, IHC-Fr, WB
免疫原	Mouse IgG, whole molecule

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.01% Sodium azide Constituents: 0.42% Potassium phosphate, 0.87% Sodium chloride
纯度	Affinity purified
纯化说明	This product was prepared from monospecific antiserum by immunoaffinity chromatography using Mouse IgG coupled to agarose beads.
克隆	多克隆
同种型	IgG

应用


The Abpromise guarantee Abpromise™承诺保证使用ab6708于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ELISA		1/150000.
IP		Use at an assay dependent concentration.
Immunodiffusion		Use at an assay dependent concentration.

应用	Ab评论	说明
Conjugation		Use at an assay dependent concentration.
ICC/IF		1/1000 - 1/5000.
IHC-P		Use at an assay dependent concentration.
IHC-Fr		Use at an assay dependent concentration.
WB		1/3000 - 1/15000.

图片



Goat Anti-Mouse IgG H&L (ab6708) at 1/1000 dilution + Mouse IgG at 0.1 µg


Secondary

Secondary antibody: DyLight™ 800 goat secondary antibody at 1/20000 dilution

Observed band size: 28,55 kDa

ab6708 incubation overnight at 4°C.

Both heavy (55kDa) and light (28kDa) chains are detected.



All lanes : Goat Anti-Mouse IgG H&L (ab6708) at 1/1000 dilution

Lane 1 : Mouse IgG (whole molecule)

Lane 2 : Mouse IgG Fc fragment

Lane 3 : Mouse IgG Fab fragment

Lane 4 : Mouse IgM kappa

Lane 5 : Mouse serum

Lysates/proteins at 0.05 µg per lane.

Secondary

All lanes : HRP Donkey Anti-Goat IgG at 1/40000 dilution

ab6708 incubated for 1 hour at room temperature.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors