

### Goat Anti-Mouse IgG H&L (HRP polymer) ab214879

★★★★★ [2 Abreviews](#) [12 References](#) [1 图像](#)

#### 概述

产品名称	山羊抗小鼠IgG H&L (HRP polymer)
宿主	Goat
靶标种属	Mouse
特异性	Based on IEP, ab214879 reacts with heavy (gamma) chains on mouse IgG, light chains on all mouse immunoglobulins. No cross-reactivity is observed to non-immunoglobulin mouse serum proteins, serum proteins from bovine, horse, human, pig or rabbit.
经测试应用	<b>适用于:</b> Immunoelectrophoresis, IHC-P
免疫原	Full length protein corresponding to Mouse IgG. (Purified Mouse IgG, whole molecule).
偶联物	HRP polymer

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: 0.1% Proclin 300 Constituents: 1% BSA, 98% PBS
纯度	ab214879 is liquid with light red color. Affinity purified
纯化说明	ab214879 was affinity purified using solid phase Mouse IgG. Purity >95% by SDS-PAGE.
克隆	多克隆
同种型	IgG

#### 应用


**The Abpromise guarantee**

**Abpromise™** 承诺保证使用ab214879于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Immunoelectrophoresis		Use at an assay dependent concentration.
IHC-P	★★★★★ (1)	Use at an assay dependent concentration. This product is provided in a pre-diluted, ready-to-use format. Incubate 150 µl per slide for 1 hour.

## 图片



Immunohistochemical analysis of Triple positive breast cancer tissue stained with Ki-67 and detected with ab214879.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Goat Anti-Mouse IgG H&L (HRP polymer) (ab214879)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors