

Goat Anti-Mouse IgG1 (HRP) preadsorbed ab98693

★★★★★ [1 Abreviews](#) [18 References](#) [1 图像](#)

概述

产品名称	山羊抗小鼠IgG1 (HRP)预吸附二抗
宿主	Goat
靶标种属	Mouse
特异性	<p>By immunoelectrophoresis and ELISA this antibody reacts specifically with mouse IgG1.</p> <p>Cross reactivity with immunoglobulins is less than 2%. No antibody was detected against non-immunoglobulin serum proteins.</p> <p>Less than 1% cross reactivity to human and rat IgG1 was detected. This antibody may cross react with IgG1 from other species.</p>
经测试应用	适用于: ICC, IHC-P, ELISA, WB
最小交叉反应	Human, Rat more details
偶联物	HRP

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 6.8 Constituents: 0.2% BSA, 0.05% CMIT/MIT based preservative
纯度	Immunogen affinity purified
纯化说明	Antiserum was solid phase adsorbed to ensure subclass specificity. Antiserum was cross adsorbed using human and rat immunosorbents to remove cross reactive antibodies. The antibody to mouse IgG1 was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to horseradish peroxidase (HRP).
克隆	多克隆
同种型	IgG

应用


Abpromise™ 承诺保证使用ab98693于以下的经测试应用

The Abpromise guarantee

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ICC		Use at an assay dependent dilution.
IHC-P		1/200 - 1/2000.
ELISA		1/10000 - 1/100000. Similar dilution to primary antibody.
WB	★★★★★ (1)	1/2000 - 1/20000. Colorimetric 1/5000 - 1/25000 Chemiluminescent

图片



IHC image of Histone H1.0 staining in human colon formalin fixed paraffin embedded tissue section*. The section was pre-treated using pressure cooker heat mediated antigen retrieval with sodium citrate buffer (pH6) for 30mins. The section was incubated with **ab11080**, 7.5µg/ml overnight at +4°C. An HRP-conjugated secondary (ab98693, 1/2000 dilution) was used for 1hr at room temperature. The section was counterstained with haematoxylin and mounted with DPX.

The inset negative control image is secondary-only at 1/500 dilution.

*Tissue obtained from the Human Research Tissue Bank, supported by the NIHR Cambridge Biomedical Research Centre

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors