

Goat Anti-Mouse IgG2a heavy chain ab9163

1 References

概述	
产品名称	山羊抗小鼠IgG2a heavy chain
宿主	Goat
靶标种属	Mouse
特异性	The affinity purified anti IgG2a only reacts with mouse IgG2a by immunoelectrophoresis and double diffusion when tested at 1 mg/ml against all mouse immunoglobulin classes and subclasses at 1 mg/ml.
经测试应用	适用于: ELISA, Immunodiffusion, Immunoelectrophoresis, WB
免疫原	Goats were immunized with highly purified mouse IgG2a in Freund's adjuvant

性能	
形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	pH: 7.20
纯度	Immunogen affinity purified
纯化说明	The antiserum was solid phase absorbed to make it IgG2a heavy chain specific and the anti IgG2A subsequently immunoaffinity purified off an IgG2a immunoabsorbant.
克隆	多克隆

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab9163于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ELISA		Use at an assay dependent concentration.
Immunodiffusion		Use at an assay dependent concentration.

应用	Ab评论	说明
Immunoelectrophoresis		Use at an assay dependent concentration.
WB		Use at an assay dependent concentration.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors