

Goat Anti-Mouse IgG2a H&L (PE) ab74490

3 References

概述

产品名称	山羊抗小鼠IgG2a H&L (PE)
宿主	Goat
靶标种属	Mouse
经测试应用	适用于: Flow Cyt
免疫原	Highly purified mouse IgG isotype in Freund's adjuvant.
偶联物	PE. Ex: 488nm, Em: 575nm

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 7.40 Preservative: 0.013% Sodium azide Constituents: 1.71% Sucrose, 0.1% BSA, 0.87% Sodium chloride, 1.64% Sodium phosphate
纯度	Size exclusion
纯化说明	Purified using size exclusion chromatography as well as some proprietary methods.
克隆	多克隆
同种型	IgG2a

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab74490于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Flow Cyt		Use at an assay dependent dilution.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors