



Product datasheet

Goat Anti-Mouse IgG2b heavy chain (HRP) preadsorbed
ab98703

3 References

概述	
产品名称	山羊抗小鼠IgG2b heavy chain (HRP)预吸附二抗
宿主	Goat
靶标种属	Mouse
经测试应用	适用于: IHC-P, ELISA, WB, ICC
最小交叉反应	Human, Rat
偶联物	HRP
more details	

性能	
形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 6.8 Constituents: 0.2% BSA, 0.05% CMIT/MIT based preservative
纯度	Immunogen affinity purified
纯化说明	Antiserum was solid phase adsorbed to ensure subclass specificity. Antiserum was cross adsorbed using human and rat immunosorbents to remove cross reactive antibodies. This antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to HRP.
克隆	多克隆
同种型	IgG

应用	
The Abpromise guarantee	Abpromise™承诺保证使用ab98703于以下的经测试应用
“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。	

应用	Ab评论	说明
IHC-P		1/200 - 1/2000.
ELISA		1/10000 - 1/100000. (Primary EISA).
WB		1/2000 - 1/25000. (Colorimetric: 1/2000-1/20000 dilution, Chemiluminescent: Use 1/5000-1/25000 dilution).
ICC		Use at an assay dependent dilution.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors