

Goat Anti-Mouse IgM (SureLight® PE) ab99917

1 图像

概述

产品名称	山羊抗小鼠IgM (SureLight® PE)
宿主	Goat
靶标种属	Mouse
经测试应用	适用于: Flow Cyt
免疫原	Full length native protein (purified) corresponding to IgM.
偶联物	SureLight® PE

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 7.20 Preservative: 0.05% Sodium azide Constituents: 0.58% Sodium chloride, 1.64% Sodium phosphate, 1.71% Sucrose, 0.1% BSA
共轭说明	Conjugated to R-Phycoerythrin. Visible absorption maxima 565>540>498 Emission maximum 578
克隆	多克隆
同种型	IgG
常规说明	Avoid exposure to heat and light.


应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab99917于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Flow Cyt		Use a concentration of 3 µg/ml.

图片



Flow Cytometry - Goat Anti-Mouse IgM (SureLight® PE) (ab99917)

Fluorescence excitation and emission spectra of R-Phycoerythrin in 100mM sodium phosphate (pH7.2) + 1mM EDTA and 1mM sodium azide. Emission scan was taken with excitation at 498nm.

Excitation scan was taken with an emission at 575nm.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors