

Goat Anti-Rabbit IgG H&L ab2892

概述

产品名称	山羊抗兔IgG H&L
宿主	Goat
靶标种属	Rabbit
特异性	This antibody cocktail reacts with heavy and light chains on rabbit IgG and mouse IgG (reacts with all IgG subtypes).
经测试应用	适用于: IHC-P, IHC-Fr, ICC

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Constituents: 0.04% Proclin 300, 0.5% Proclin 950, 99% Water
共轭说明	HRP or AP conjugated to a polymer backbone, which is collectively conjugated to the antibody. The end result is two antibody subsets: 1) Multiple HRP molecules linked to goat anti-rabbit antibodies. 2) Multiple AP molecules linked to goat anti-mouse antibodies.
克隆	多克隆
同种型	IgG
常规说明	<p>This product is specially designed for a rapid double stain procedure. The overall staining procedure can be done in 4 steps and in less than two hours. Avidin-biotin blocking procedures are not necessary when using this product. This new and innovative AP and HRP polymerization technology provides a significant increase in staining sensitivity when compared to other conventional AP or HRP conjugated secondary antibodies.</p> <p>This product is provided in a pre-diluted, ready-to-use format. Procedure (adapt in-house protocols to accommodate): 1) A primary antibody cocktail with a rabbit monoclonal/polyclonal and a mouse monoclonal/polyclonal antibody must be applied for the standard duration. Rinse in TBS. (DO NOT USE PBS). 2) Apply ab2892 for 20 - 30 minutes then rinse in TBS (DO NOT USE PBS). 3) Chromogen 1: Develop AP using suitable chromogen. Rinse in TBS. 4) Chromogen 2: Develop HRP using suitable chromogen. Rinse in TBS. Equally suited to both manual and automated immunostaining procedures. This product has not been tested in any application requiring further dilution of the antibody and therefore it is not recommended for use in such applications.</p>

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab2892于以下的经测试应用

“应用说明”部分下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IHC-P		1/1.
IHC-Fr		1/1.
ICC		1/1.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors