

Goat Anti-Rabbit IgG H&L (FITC) preadsorbed ab7086

★★★★★ 1 Abreviews 23 References

概述	
产品名称	山羊抗兔IgG H&L (FITC)预吸附二抗
宿主	Goat
靶标种属	Rabbit
经测试应用	适用于: ICC/IF, Immunomicroscopy, Flow Cyt, IHC-P, IHC-Fr, ELISA
最小交叉反应	Chicken, Cow, Goat, Guinea pig, Hamster, Horse, Human, Mouse, Rat, Sheep <a href="#">more details</a>
免疫原	Full length native Rabbit IgG (purified).
偶联物	FITC. Ex: 493nm, Em: 528nm
性能	
形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 6.50 Preservative: 0.01% Sodium azide Constituents: 0.42% Tripotassium orthophosphate, 0.87% Sodium chloride, 1% BSA
纯度	Affinity purified
纯化说明	This product was prepared from monospecific antiserum by immunoaffinity chromatography using Rabbit IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities.
共轭说明	Fluorescein isothiocyanate (FITC) (Molecular Weight 390 daltons) Absorption Wavelength: 495 nm Emission Wavelength: 528 nm Fluorochrome/Protein Ratio: 2.7 moles FITC per mole of Goat IgG
克隆	多克隆
同种型	IgG
应用	

The Abpromise guarantee Abpromise™承诺保证使用ab7086于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ICC/IF		
Immunomicroscopy		
Flow Cyt		
IHC-P		
IHC-Fr	★★★★★ (1)	
ELISA		

#### 应用说明

Flow Cyt: Use at an assay dependant concentration.

IHC-P: Use at an assay dependant concentration.

IHC-Fr: Use at an assay dependant concentration.

IM: Use at an assay dependant concentration.

Not tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

#### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

#### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors