

# Goat Anti-Rat IgG H&L (Cy3<sup>®</sup>) preadsorbed ab6953

## 17 References

### 概述

产品名称	山羊抗大鼠IgG H&L (Cy3 <sup>®</sup> )预吸附二抗
宿主	Goat
靶标种属	Rat
经测试应用	适用于: ICC/IF, ICC, Flow Cyt, IHC-P, IHC-Fr, ELISA
最小交叉反应	Chicken, Cow, Goat, Guinea pig, Hamster, Horse, Human, Mouse, Rabbit, Sheep <a href="#">more details</a>
免疫原	Full length native Rat IgG (purified).
偶联物	Cy3 <sup>®</sup> . Ex: 552nm, Em: 565nm

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 6.50 Preservative: 0.01% Sodium azide Constituents: 0.42% Tripotassium orthophosphate, 0.87% Sodium chloride, 1% BSA
纯度	Immunogen affinity purified
纯化说明	This product was prepared from monospecific antiserum by immunoaffinity chromatography using Rat IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities.
共轭说明	Cy3.29 (Cyanine 3.29-OSu) (Molecular Weight 949 daltons) Absorption Wavelength: 552 nm Emission Wavelength: 565 nm Fluorochrome/Protein Ratio: 8.5 moles Cy3 per mole of Goat IgG
克隆	多克隆
同种型	IgG
常规说明	Cy <sup>™</sup> and CyDye <sup>™</sup> are registered trademarks of Cytiva.  This secondary antibody is specifically designed for the detection of multiple primary antibodies (polyclonal or monoclonal) of different host species in experiments where cells are simultaneously labeled without unwanted cross reaction.

## 应用

---

### The Abpromise guarantee **Abpromise™**承诺保证使用ab6953于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ICC/IF		1/1000 - 1/5000.
ICC		Use at an assay dependent concentration.
Flow Cyt		Use at an assay dependent dilution.
IHC-P		Use at an assay dependent concentration.
IHC-Fr		Use at an assay dependent concentration.
ELISA		1/10000 - 1/50000.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors