

Goat Anti-Rat IgG H&L (HRP) preadsorbed ab7097

54 References 2 图像

概述	
产品名称	山羊抗大鼠IgG H&L (HRP)预吸附二抗
宿主	Goat
靶标种属	Rat
经测试应用	适用于: ICC/IF, IHC-P, IHC-Fr, WB, ELISA, Dot blot
最小交叉反应	Chicken, Cow, Goat, Guinea pig, Hamster, Horse, Human, Mouse, Rabbit, Sheep <a href="#">more details</a>
免疫原	Full length native Rat IgG (purified).
偶联物	HRP

性能	
形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 6.50 Preservative: 0.01% Gentamicin sulphate Constituents: 0.42% Tripotassium orthophosphate, 0.87% Sodium chloride, 1% BSA
纯度	Affinity purified
纯化说明	This product was prepared from monospecific antiserum by immunoaffinity chromatography using Rat IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities.
共轭说明	Horseradish Peroxidase (HRP)
克隆	多克隆
同种型	IgG


应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab7097于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ICC/IF		1/500 - 1/1500.
IHC-P		Use at an assay dependent dilution.
IHC-Fr		Use at an assay dependent dilution.
WB		Use at an assay dependent dilution.
ELISA		1/50000.
Dot blot		Use at an assay dependent dilution.

## 图片



Goat Anti-Rat IgG H&L (HRP) preadsorbed (ab7097) at 1/1000 dilution + Rat IgG, 50 ng

Both the heavy and light chains of rat IgG are detected.



Whole mount immunostaining of E11.5 mouse embryo (Blood vessels in the embryonic head). This staining was performed with rat monoclonal antibody to Pecam followed by staining with ab7097 (HRP conjugated).

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors