

### Goat Anti-Rat IgM mu chain (Biotin) ab97178

#### 4 References

#### 概述

产品名称	山羊抗大鼠IgM mu chain (Biotin)
宿主	Goat
靶标种属	Rat
经测试应用	适用于: IHC-P, ICC/IF, ELISA, WB
偶联物	Biotin

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	pH: 6.8 Preservative: 0.09% Sodium azide Constituents: 0.2% BSA, PBS
纯度	Immunogen affinity purified
纯化说明	Antiserum was solid phase adsorbed to ensure class specificity. This antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to Biotin. Biotinylated antibody was demonstrated with avidin/oxidase when coated on microtiter wells.
共轭说明	Biotin/ antibody protein ratio is 4:1
克隆	多克隆
同种型	IgG

#### 应用

**The Abpromise guarantee**      **Abpromise™** 承诺保证使用ab97178于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IHC-P		1/250 - 1/2500.

应用	Ab评论	说明
ICC/IF		1/100 - 1/500.
ELISA		1/10000 - 1/200000. (Primary).
WB		1/10000 - 1/200000. (Chemiluminescent).

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors