



Product datasheet

Anti-Human IgG4 pFc' [HP6023] (PE) ab99819

概述

产品名称	抗人IgG4 pFc' [HP6023] (PE)
宿主	Mouse
靶标种属	Human
特异性	ab99819 reacts with the pFc' portion of the heavy chain of IgG4.
经测试应用	适用于: ELISA, Dot blot, IHC-Fr, ICC/IF, WB
免疫原	Full length protein corresponding to IgG4. Human IgG4 myeloma protein
偶联物	PE. Ex: 488nm, Em: 575nm

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C. Store In the Dark.
存储溶液	pH: 7.30 Preservative: 0.093% Sodium azide Constituents: 98% PBS, 1% Sucrose
纯度	Protein A purified
克隆	单克隆
克隆编号	HP6023
同种型	IgG3
轻链类型	kappa

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab99819于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ELISA		Use at an assay dependent dilution.

应用	Ab评论	说明
Dot blot		Use at an assay dependent dilution. for Dot and slot immunoblotting
IHC-Fr		Use at an assay dependent dilution.
ICC/IF		Use 0.1µg for 10 ⁶ cells.
WB		Use at an assay dependent dilution.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors