

Mouse monoclonal [hlgA1] Anti-Human IgA H&L ab433

1 References

概述	
产品名称	小鼠单克隆抗体[hlgA1]抗人IgA H&L
宿主	Mouse
靶标种属	Human
特异性	hlgA1 specifically recognizes human IgA. hlgA - 100% hlgA (secretory) - 107% hlgA (myeloma) - 135% hlgA1 (myeloma) - 145% hlgA2 (myeloma) - 91% hlgA1 (myeloma) - 104% hlgA2 (myeloma) - 200% hlgA alpha-chain - 138% hlgG - <0.3% hlgG - <0.3% hlgM (plasma) - 2.6% hlgM (Waldenström's macroglobulinaemia) - <0.3% Affinity constant: 3x10 ⁹ L/mol human plasma IgA.
经测试应用	适用于: ELISA, RIA
性能	
形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	pH: 6.00 Preservative: 0.095% Sodium azide Constituents: 0.335% PBS, 1.0878% Sodium citrate, 0.9% Sodium chloride
纯度	Protein A purified
克隆	单克隆
克隆编号	hlgA1
同种型	IgG1
轻链类型	unknown
应用	

The Abpromise guarantee **Abpromise™**承诺保证使用ab433于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ELISA		Use at an assay dependent concentration.
RIA		Use at an assay dependent concentration.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors