

Mouse monoclonal [MARM-4] Anti-Rat IgM H&L (HRP) ab11678

概述

产品名称	小鼠单克隆抗体[MARM-4]抗大鼠IgM H&L (HRP)
宿主	Mouse
靶标种属	Rat
特异性	Ab11678 recognises rat IgM and does not react with other rat immunoglobulin isotypes.
经测试应用	适用于: ELISA
免疫原	Full length protein corresponding to Rat IgM.
偶联物	HRP

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: 0.09% Sodium azide Constituents: PBS, 50% Glycerol
纯化说明	Purified IgG prepared by affinity chromatography.
克隆	单克隆
克隆编号	MARM-4
骨髓瘤	Sp2/0
同种型	IgG1

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab11678于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ELISA		

应用说明 ELISA: Use at a concentration of 0.25 µg/ml as detecting antibody.

Not tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors