

Mouse monoclonal [SB62a] Anti-Rabbit IgG light chain ab99695

概述

| | |
|-------|---------------------------------|
| 产品名称 | 小鼠单克隆抗体[SB62a]抗兔IgG light chain |
| 宿主 | Mouse |
| 靶标种属 | Rabbit |
| 经测试应用 | 适用于: WB, ELISA |
| 免疫原 | Rabbit Light Chain |

性能

| | |
|------|---|
| 形式 | Liquid |
| 存放说明 | Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle. |
| 存储溶液 | pH: 8.20 Constituent: 100% Borate buffered saline |
| 纯度 | Protein A purified |
| 克隆 | 单克隆 |
| 克隆编号 | SB62a |
| 同种型 | IgG2b |
| 轻链类型 | lambda |

应用

The Abpromise guarantee **Abpromise™** 承诺保证使用ab99695于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

| 应用 | Ab评论 | 说明 |
|-------|------|-------------------------------------|
| WB | | Use at an assay dependent dilution. |
| ELISA | | Use a concentration of 1 µg/ml. |

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors