

Rabbit Anti-Goat IgG H&L (beta-galactosidase) ab136712

3 References

概述	
产品名称	兔抗山羊IgG H&L (beta-galactosidase)
宿主	Rabbit
靶标种属	Goat
经测试应用	适用于: IHC-Fr, IHC-P, ICC, WB, ELISA
免疫原	Highly purified whole immunoglobulin from Goat.
偶联物	beta-galactosidase

性能	
形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Constituents: 49% PBS, 50% Glycerol
纯度	Immunogen affinity purified
纯化说明	ab136712 is purified by a very stringent purification protocol where non-specific and specific purifications alternate (affinity chromatography).
克隆	多克隆
同种型	IgG
常规说明	Preparation of Substrate Solutions: 0.1M phosphate buffer pH7, 1mM MgSO ₄ , 0.2mM MnSO ₄ , 2mM Mg titriplex. Add extemporaneously 7µl/ml beta-mercaptoethanol and 0.9mg/ml ortho-nitrophenyl-beta-galactopyranoside. Incubate for 30 minutes at 37°C. Stop the reaction with V/2 volume of 1M Na ₂ CO ₃ .

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab136712于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明

应用	Ab评论	说明
IHC-Fr		Use at an assay dependent concentration.
IHC-P		Use at an assay dependent concentration.
ICC		Use at an assay dependent concentration.
WB		Use at an assay dependent concentration.
ELISA		1/100 - 1/1000.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors