

Rabbit Anti-Human IgG H&L (Alkaline Phosphatase) ab6760

4 References

概述	
产品名称	兔抗人IgG H&L (Alkaline Phosphatase)
宿主	Rabbit
靶标种属	Human
经测试应用	适用于: Immunomicroscopy, WB, ICC/IF, IHC-Fr, Dot blot, IHC-P, ELISA
免疫原	human IgG whole molecule
偶联物	Alkaline Phosphatase
性能	
形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: 0.1% Sodium azide Constituents: 0.878% Sodium chloride, 1% BSA, 0.788% Tris HCl, 0.0014% Zinc chloride, 0.0095% Magnesium chloride, 50% Glycerol
纯度	Affinity purified
纯化说明	This product was prepared from monospecific antiserum by immunoaffinity chromatography using human IgG coupled to agarose beads.
共轭说明	Alkaline Phosphatase (Calf Intestine) (Molecular Weight 140,000 daltons)
克隆	多克隆
应用	

The Abpromise guarantee      **Abpromise™**承诺保证使用ab6760于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Immunomicroscopy		
WB		

应用	Ab评论	说明
ICC/IF		
IHC-Fr		
Dot blot		
IHC-P		
ELISA		

**应用说明** suitable for immunoblotting (western or dot blot), ELISA and immunohistochemistry as well as other phosphatase-antibody based enzymatic assays requiring lot-to-lot consistency

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors