

# Rabbit Anti-Human IgG H&L (FITC) ab6755

## 2 References

### 概述

|       |   |
|-------|---|
| 产品名称  | 兔抗人IgG H&L (FITC)   |
| 宿主    | Rabbit  |
| 靶标种属  | Human   |
| 经测试应用 | <b>适用于:</b> Immunomicroscopy, Flow Cyt, ICC/IF, IHC-Fr, ELISA |
| 免疫原   | human whole molecule  |
| 偶联物   | FITC. Ex: 493nm, Em: 528nm                                    |

### 性能

|      |  |
|------|--|
| 形式   | Liquid   |
| 存放说明 | Shipped at 4°C. Store at +4°C.   |
| 存储溶液 | pH: 6.50<br>Preservative: 0.01% Sodium azide<br>Constituents: 0.878% Sodium chloride, 1% BSA, 0.424% Tripotassium orthophosphate   |
| 纯度   | Affinity purified  |
| 纯化说明 | This product was prepared from monospecific antiserum by immunoaffinity chromatography using human IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities. |
| 共轭说明 | Diaminotriazinylaminofluorescein (DTAF) (Molecular Weight 530 daltons) Absorption Wavelength: 495 nm Emission Wavelength: 528 nm Fluorochrome/Protein Ratio: 2.3 moles DTAF per mole of Rabbit IgG         |
| 克隆   | 多克隆  |

### 应用

**The Abpromise guarantee**      **Abpromise™** 承诺保证使用ab6755于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

| 应用               | Ab评论 | 说明                                       |
|------------------|------|--|
| Immunomicroscopy |      | 1/1000 - 1/5000.                         |
| Flow Cyt         |      | 1/500 - 1/2500.                          |
| ICC/IF           |      | 1/1000 - 1/5000.                         |
| IHC-Fr           |      | Use at an assay dependent concentration. |
| ELISA            |      | 1/10000 - 1/50000.                       |

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors