

Rabbit Anti-Mouse IgG+IgM+IgA H&L (FITC) ab8517

12 References 1 图像


概述	
产品名称	兔抗小鼠IgG+IgM+IgA H&L (FITC)
宿主	Rabbit
靶标种属	Mouse
特异性	Antibodies were raised in rabbit to mouse IgG, IgM and IgA
经测试应用	适用于: ICC/IF, IHC-P
偶联物	FITC. Ex: 493nm, Em: 528nm
性能	
形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 7.50 Constituents: 0.164% Sodium phosphate, 50% Glycerol (glycerin, glycerine), 0.58% Sodium chloride
纯度	Immunogen affinity purified
纯化说明	Specific antibodies were absorbed by incubation with Sepharose-bound mouse immunoglobulins. Specific antibodies were eluted by acidic buffer at pH 2.5 followed by neutralisation and dialysis. After repeated binding with immobilized mouse immunoglobulins a minimum of 65% protein bound.
共轭说明	The purified polyclonal was conjugated to fluorescein isothiocyanate isomer I followed by gel-filtration and ion-exchange chromatography to clear unbound conjugate.
克隆	多克隆
常规说明	F/P ratio is 3-5 mol Fluorescein per 1 mol of rabbit Ig. There is low background fluorescence of normal tissue in contrast to specific staining of sections treated with monoclonal mouse anti-human collagen antibodies
应用	

The Abpromise guarantee      Abpromise™承诺保证使用ab8517于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ICC/IF		Use at an assay dependent dilution. PubMed: 15483223
IHC-P		Use at an assay dependent dilution.

## 图片



Ab8517 was used at dilution 1/200 with the primary antibody **ab18230** in ICC/IF. See the review on **ab18230**.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors