

Rabbit Anti-Sheep IgG H&L (Alkaline Phosphatase) preadsorbed ab7112

概述

产品名称	兔抗羊IgG H&L (Alkaline Phosphatase)预吸附二抗
宿主	Rabbit
靶标种属	Sheep
经测试应用	适用于: ICC/IF, IHC-P, IHC-Fr, WB, ELISA, Immunomicroscopy, Dot blot
最小交叉反应	Human
免疫原	Full length native Sheep IgG (purified).
偶联物	Alkaline Phosphatase

[more details](#)

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.01% Sodium azide Constituents: 0.00136% Zinc chloride, 0.0095% Magnesium chloride, 0.79% Tris HCl, 50% Glycerol (glycerin, glycerine), 0.87% Sodium chloride, 1% BSA
纯度	Affinity purified
纯化说明	Alkaline Phosphatase (Calf Intestine) (Molecular Weight 140,000 daltons)
共轭说明	This product was prepared from monospecific antiserum by immunoaffinity chromatography using Sheep IgG coupled to agarose beads
克隆	多克隆

应用

The Abpromise guarantee Abpromise™承诺保证使用ab7112于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ICC/IF		1/200 - 1/1000.
IHC-P		Use at an assay dependent dilution.
IHC-Fr		Use at an assay dependent dilution.
WB		1/400 - 1/1600.
ELISA		1/1000 - 1/5000.
Immunomicroscopy		Use at an assay dependent dilution.
Dot blot		Use at an assay dependent dilution.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors