

Rabbit Anti-Sheep IgG H&L preadsorbed ab7106

概述

产品名称	兔抗羊IgG H&L预吸附二抗
宿主	Rabbit
靶标种属	Sheep
经测试应用	适用于: WB, ELISA, ICC/IF, IHC-P, IP, Immunodiffusion, Conjugation, IHC-Fr
最小交叉反应	Human
免疫原	Full length native Sheep IgG (purified).

[more details](#)

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.01% Sodium azide Constituents: 0.164% Sodium phosphate, 1.45% Sodium chloride
纯度	Affinity purified
纯化说明	This product was prepared from monospecific antiserum by immunoaffinity chromatography using Sheep IgG coupled to agarose.
克隆	多克隆

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab7106于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		
ELISA		
ICC/IF		

应用	Ab评论	说明
IHC-P		
IP		
Immunodiffusion		
Conjugation		
IHC-Fr		

应用说明

Conjugation: Use at an assay dependant concentration.

ID: Use at an assay dependant concentration.

IHC-P: Use at an assay dependant concentration.

IHC-Fr: Use at an assay dependant concentration.

IP: Use at an assay dependant concentration.

Not tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors