

Rabbit Anti-Syrian Hamster IgG H&L (HRP) ab6783

1 References

概述

产品名称	兔抗叙利亚仓鼠IgG H&L (HRP)
宿主	Rabbit
靶标种属	Syrian hamster
经测试应用	适用于: ICC, Dot blot, IHC-P, IHC-Fr, Immunomicroscopy, WB, ICC/IF, ELISA
免疫原	hamster IgG whole molecule
偶联物	HRP

性能

形式	Lyophilized: Store vial at 4° C prior to restoration. Restore with 1.0 mL of deionized water (or equivalent). For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. Expiration date is one (1) year from date of restoration.
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: 0.01% Gentamicin sulphate Constituents: 0.87% Sodium chloride, 1% BSA, 0.42% Potassium phosphate
纯度	IgG fraction
纯化说明	Horseradish Peroxidase (HRP).
共轭说明	This product was prepared from monospecific antiserum by immunoaffinity chromatography using hamster IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities
克隆	多克隆

应用

The Abpromise guarantee

Abpromise™ 承诺保证使用 ab6783 于以下的经测试应用

“应用说明”部分下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ICC		Use at an assay dependent dilution.
Dot blot		Use at an assay dependent concentration.
IHC-P		Use at an assay dependent concentration.
IHC-Fr		Use at an assay dependent concentration.
Immunomicroscopy		Use at an assay dependent concentration.
WB		1/1000 - 1/5000.
ICC/IF		1/500 - 1/2500.
ELISA		1/16000 - 1/64000.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors