

## Product datasheet

# Rat F(ab')<sub>2</sub> monoclonal [SB73a] Anti-Mouse IgM mu chain ab99600

### 概述

产品名称	大鼠F(ab') <sub>2</sub> 单克隆抗体[SB73a]抗小鼠IgM mu chain
宿主	Rat
靶标种属	Mouse
特异性	ab99600 reacts with the $\mu$ chain of BALB/c Mouse IgM as demonstrated by ELISA.
经测试应用	适用于: Flow Cyt, ELISA, ICC/IF

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
存储溶液	pH: 8.20 Constituent: 100% Borate buffered saline  No preservatives or amine-containing buffer salts added.
纯度	Protein A purified
克隆	单克隆
克隆编号	SB73a
同种型	IgG1

### 应用

**The Abpromise guarantee**      **Abpromise™** 承诺保证使用ab99600于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Flow Cyt		Use 1 $\mu$ g for 10 <sup>6</sup> cells.
ELISA		Use a concentration of 1 $\mu$ g/ml.

应用	Ab评论	说明
ICC/IF		Use at an assay dependent dilution.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors