

Product datasheet

Anti-v-Myb + a-Myb antibody [5E11] ab10934

1 References

概述

产品名称	Anti-v-Myb + a-Myb抗体[5E11]
描述	小鼠单克隆抗体[5E11] to v-Myb + a-Myb
宿主	Mouse
特异性	ab10934 recognises v-myb and A-myb from Xenopus cells. The antibody does not recognise c-myb in human or mouse.
经测试应用	适用于: IHC-Fr, WB, IP
种属反应性	与反应: Xenopus laevis 不与反应: Mouse, Human
免疫原	Recombinant v-myb protein, DNA binding protein.
阳性对照	Xenopus testes lysate or tissue sections.
常规说明	This product should be stored undiluted. Storage in frost free freezers is not recommended. Should this product contain a precipitate we recommend microcentrifugation before use.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: 0.09% Sodium Azide Constituents: PBS
纯度	Protein G purified
纯化说明	Purified from tissue culture supernatant.
克隆	单克隆
克隆编号	5E11
骨髓瘤	Sp2/0
同种型	IgG2a

应用

Our [Abpromise guarantee](#) covers the use of **ab10934** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
IHC-Fr		Use at an assay dependent concentration.
WB		Use at an assay dependent concentration. Predicted molecular weight: 72 kDa.
IP		Use at an assay dependent concentration.

靶标

相关性 Transcriptional activator; DNA-binding protein that specifically recognize the sequence 5'-YAAC[GT]G-3'. Plays an important role in the control of proliferation and differentiation of hematopoietic progenitor cells.

细胞定位 Nuclear

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors