

Product datasheet

Anti-whole Human Serum antibody (HRP) ab34671

概述

产品名称	Anti-whole人Serum抗体(HRP)
描述	兔多克隆抗体to whole人Serum (HRP)
宿主	Rabbit
偶联物	HRP
经测试应用	适用于: WB, ELISA, Dot blot, Electron Microscopy
种属反应性	与反应: Human
免疫原	Human serum proteins.

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: None Constituents: 10mg/ml BSA (IgG and Protease free), 0.14M Sodium chloride, 0.01M Sodium phosphate. pH 7.4
纯度	IgG fraction
纯化说明	This product is purified by a multi-step process including delipidation, salt fractionation and ion exchange chromatography followed by extensive dialysis against buffer.
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab34671** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		
ELISA		

应用	Ab评论	说明
Dot blot		
Electron Microscopy		

应用说明

Dot: Use at an assay dependent dilution.
 Immunoperoxidase Electron Microscopy: Use at an assay dependent dilution.
 ELISA: Use at an assay dependent dilution.
 WB: Use at an assay dependent dilution.
 Suitable for other peroxidase-antibody based assays.

Not yet tested in other applications.
 Optimal dilutions/concentrations should be determined by the end user.

靶标

相关性

The human serum proteome as a whole reflects changes in the physiological state of the human body. The onset of diseases invariably induces changes in the protein content, perhaps by producing more of one protein, less of another, or by producing a whole new protein. In all of these situations, the changes offer a direct route for diagnosing the disease.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors